



WORK INSTRUCTIONS FOR SOLVENTS

Tema:
Arbeidsinstruks
Revised date:
18.10.23
Side:
1 av 1

H220	Extremely flammable gas
H221	Flammable gas
H224	Extremely flammable liquid and vapour
H225	Highly flammable liquid and vapour
H242	Heating may cause a fire

PURCHASING

- The person in charge of purchase on the unit or in the research group

UTILIZATION

- Working with chemicals labels flammable, highly flammable and extremely flammable should take place in a fume cupboard.
- Read the material safety data sheet before starting to work with new chemicals.
- During the first five months of the pregnancy, pregnant women should not work with this kind of chemicals. A separate work instruction for pregnant and breast feeding should be conferred.

PROTECTIVE EQUIPMENT

- Use laboratory coat, protective shoes, splash goggles and protective gloves when working with corrosive chemicals. For most corrosive chemicals Sol-Vex gloves is the correct glove. You can find a summary of the most common chemicals and preferred glove [here](#).

DISPOSAL

- Follow [UIB's procedure for disposal](#).
- Use the blue jar for disposal of solvents with halogens. Disposal group 7151.
- Use the transparent jar for disposal of solvents without halogens. Disposal group 7152.
- All departments/groups also have their own disposal contact.

LOCAL PROCEDURES

- The procedure must be adapted to the individual unit.